

| Features & Benefits | |
|-------------------------|---|
| Basic Navigation | High Precision Electromagnetic tracking with sensor at instrument tip |
| | Malleable Instruments Manually adjustment for patient anatomy |
| | Plug & Play No recalibration of Scopis instruments necessary (except universal tracker) |
| | OR Integration Small console easy to integrate |
| | Intuitive user-interface |
| | Full-HD Documentation DICOM import from PACS |

| Category | Product | Art. Code | Comment |
|-------------------------|--|--------------|--|
| Capital System | Electromagnetic (EM) Navigation Unit | 8000-010-003 | Mandatory |
| | Field Generator Mounting Arm | 8000-010-005 | |
| | Field Generator | 8000-010-004 | |
| Patient Tracker | Patient tracker tabs | 8000-100-001 | Mandatory |
| | Patient Tracker EM (unlimited Use) | 8000-040-001 | One of both mandatory |
| | Patient Tracker EM - 10 uses | 8000-040-002 | |
| Pointer | Precision Pointer EM, 1,5mm Ø | 8000-050-001 | One of both mandatory |
| | Pointer EM, 2,7mm Ø | 8000-050-002 | |
| Cables | S-Video cable 2.0 m | 8000-010-040 | 2 cables necessary - 1 for monitor - 1 for endo tower - Usually both DVI - Check with Customer |
| | DVI-D cable 2.0 m | 8000-010-042 | |
| | DVI to VGA cable 2.0 m | 8000-010-044 | |
| | Coaxial cable 5.0 m | 8000-010-045 | |
| | DB9 breakout cable | 8000-010-048 | |
| | DVI-I to HDMI + Component breakout cable | 8000-010-049 | |
| Mouse | Mouse, wired | 8000-030-020 | One of both mandatory |
| | Mouse, wireless | 8000-030-021 | |
| Keyboard | Medical keyboard U.S./international | 8000-030-010 | Mandatory (choose language) |
| | Medical keyboard GER | 8000-030-011 | |
| | Medical keyboard UK | 8000-030-012 | |
| | Medical keyboard ES | 8000-030-013 | |
| | Medical keyboard FR | 8000-030-014 | |
| Medical keyboard Nordic | 8000-030-015 | | |
| Software | Scopis ENT Software | 8000-020-001 | Mandatory |

| | |
|--------------------------------|--|
| Advanced (TGS & AR) | Planning Software "Building Blocks" Draw blocks according to IFAC, draw pathway |
| | Augmented Reality Guidance to target with visual overlay of planned pathway |
| | Endoscope Tracking Endoscope tracker and adapter system for different endoscopes |

| | | | |
|---|------------------------------|--------------|--|
| Endoscope Tracking | Calibration Body EM | 8000-060-002 | Mandatory |
| | Endoscope Tracker EM* | 8000-060-001 | End. Tracker* |
| | Universal Tracker EM | 8000-060-006 | or |
| | Instrument Clamp, 2-6 mm | 8000-060-011 | Universal Tracker + Clamp |
| <i>*only for 4mm x 180mm Storz Endoscopes</i> | | | |
| Software | Scopis ENT Software with TGS | 8000-020-002 | Mandatory (replaces Standard Software) |

| | |
|---------------------------|---|
| Optional Equipment | Non sterile patient registration |
| | Suction Tubes Malleable suction tubes for precise navigation |
| | Instrument Tracking Conventional instruments, forceps, shaver, suction etc.) |
| | Stand Alone System If navigation system is not integrated in Endo Tower or Operating Room |
| | Extension cable for more flexibility |

| | | | |
|---------------------|----------------------------------|--------------|---|
| Pointer | Registration Pointer EM, 2,5mm Ø | 8000-050-003 | Optional |
| Suction | Suction Tube Frazier EM | 8000-050-005 | Optional |
| | Suction Tube Eicken EM | 8000-050-006 | |
| Instrument Tracking | Universal Tracker EM | 8000-060-006 | Optional (Universal Tracker necessary for clamp usage) |
| | Instrument Clamp, Forceps | 8000-060-010 | |
| | Instrument Clamp, 2-6 mm | 8000-060-011 | |
| | Instrument Clamp, 6-10 mm | 8000-060-012 | |
| Cart | Equipment Cart PRO | 8000-030-002 | Optional |
| | HD Monitor VisionPro | 240-031-020 | Optional |
| | 4K Monitor | 240-031-050 | (one of both) |
| Cable | Navigation Tool Extension cable | 8000-050-011 | Optional |